



Manufacturer of furnaces and equipment for laboratories as well as the metal and glass industry



**For laboratories > Tubular type RSD**

**SPLIT TUBULAR FURNACE ON STAND**



Type PID	<b>RSD 60x780/25</b>
Max. temp. [°C]	<b>250</b>
Reactor diameter [mm]	<b>65</b>
Heating zone length [mm]	<b>780</b>
Supply voltage [v]	<b>230</b>
Heating zones	<b>4</b>
Thermocouple type	<b>K</b>
Power installed [kW]	<b>1,6</b>

If you have any questions please contact us:

**tel. 32 47 07 495**

**tel./fax 32 47 07 502**

**e-mail: [czylok@czylok.com.pl](mailto:czylok@czylok.com.pl)**